

Pressure Tube – *with plunger valve*



Heavy wall tube with bushing and plunger valve that allows purging of the tube. Closed bottom plunger has a hole inside, that when positioned in relation to O-Ring seal, will open Tube to the atmosphere (i.e.,

pull to close, push to open). Plunger has step-down extension that acts as stop in bushing. Complete item consists of Tube (#7 or #15 Ace-Thred), PTFE Bushing with FETFE® O-Ring, and Plunger Valve.

Length (below thread), cm	Body O.D., mm	Approx. Capacity, mL	Tube only Order Code	Complete Order Code
#7 Ace-Thred				
10.2 (4")	8	1	8648-40	8648-60
10.2 (4")	13	4	8648-41	8648-61
10.2 (4")	19	9	8648-42	8648-62
20.3 (8")	8	2	8648-45	8648-64
20.3 (8")	13	8	8648-46	8648-65
20.3 (8")	19	18	8648-47	8648-66
30.5 (12")	8	3	8648-50	8648-68
30.5 (12")	13	15	8648-51	8648-69
30.5 (12")	19	27	8648-52	8648-70
#15 Ace-Thred				
10.2 (4")	25.4	15	8648-23	8648-75
10.2 (4")	38.1	60	8648-24	8648-76
17.8 (7")	25.4	35	8648-26	8648-78
17.8 (7")	38.1	100	8648-27	8648-79
20.3 (8")	25.4	38	8648-29	8648-82
20.3 (8")	38.1	120	8648-30	8648-83
30.5 (12")	25.4	60	8648-32	8648-85
30.5 (12")	38.1	185	8648-33	8648-86
Replacement O-Rings				
Size -008 FETFE, 12/pk			7855-704	
Size -013 FETFE, 12/pk			7855-710	

Pressure Tube – *with plunger valve and thermowell*



Rugged, heavy wall tube, similar to listing above, except glass plunger valve has a 5mm I.D. thermowell for a thermocouple wire or sheathed probe up to 1/8" (3.2mm) O.D. Closed bottom plunger has a hole inside that when positioned in relation to O-Ring seal, allows Tube to be open or closed to the atmosphere (i.e., pull to close, push to

open). Step-down extension on upper end of plunger acts as stop in bushing. Complete item consists of Tube (#15 Ace-Thred), PTFE Bushing with FETFE® O-Ring, and Plunger Valve with Thermowell. See listing above for replacement Tube only with #15 Ace-Thred and 8648-19 Bushing.

Length (below thread), cm	Body O.D., mm	Approx. Capacity, mL	Tube only Order Code	Complete Order Code
#15 Ace-Thred				
10.2 (4")	25.4	15	8648-23	8648-164
10.2 (4")	38.1	48	8648-24	8648-167
17.8 (7")	25.4	35	8648-26	8648-170
17.8 (7")	38.1	100	8648-27	8648-172
20.3 (8")	25.4	38	8648-29	8648-175
20.3 (8")	38.1	120	8648-30	8648-177
30.5 (12")	25.4	60	8648-32	8648-179
30.5 (12")	38.1	185	8648-33	8648-181
Replacement O-Rings				
Size -013 FETFE, 12/pk			7855-710	